

Laser Spectroscopy: Measuring Materials' Optical Properties

Shreya Sutariya (shreyas@uchicago.edu), ERC 568

We use the **Frequency-tunable Laser Source (FLS)** to study the optical properties of materials such as absorbers, filters, anti-reflective coatings, etc that are used inside CMB telescopes. Understanding these materials as accurately as possible help us know which optical systematics can arise. We can characterize materials by taking reflection and transmission measurements on the optics bench and fitting to models to get parameters such as the loss tangent, index of refraction, etc.

Today's lab will center on reflection measurements of a common plastic called high-density polyethelene (HDPE). After measuring and plotting the measurement, we will plot its simulated performance and compare. Finally, we fit to get the exact index of refraction n of the material. (Time permitting we will also measure: a Low-Pass Edge (LPE) filter and non-ARC alumina).

Today's Objectives:

- Background: Understand measurement instrument and reflection setup.
- The Exercises are in `Notebook2_Reflection.ipynb`
- *Exercise 1:* Take reflection measurements of our samples, note their thicknesses, and the angle of incidence in the setup.
- *Exercise 2:* Simulate what reflection of such a material should look like.
- *Exercise 3:* Fit our measurement to the model to get the index of refraction of our sample.
 - Bonus: We can understand a material more precisely by knowing its complex index of refraction, $\tilde{n} = n + ik$, where the imaginary component k determines the absorptivity of the material (and the real part n is what we normally refer to as the index of refraction). For code that can fit to both the real and imaginary parts, use the FLS Lab Computer to do additional fits using `Notebook3_Bonus_complex_index_fits.ipynb`.

Background: The Frequency-tunable Laser Source (FLS)

The FLS is a broadband frequency source that is capable of emitting from ~ 20 GHz to 1200 GHz. This photomixing source and receiver package operates by using two infrared diode lasers, which emit at slightly different frequencies: ν_1 and ν_2 . The difference between the two, $\Delta\nu = \nu_2 - \nu_1$, is our desired frequency, also known as the beat frequency. The two laser frequencies are incident on a photomixer and induce an oscillating photocurrent at our beat frequency in the semiconductor material of the photomixer. A log-spiral antenna then emits this photocurrent as a THz wave.

Just as the laser system uses a photomixer to transmit a THz wave, it also uses a photomixer to detect it. This process involves an induced photocurrent in the receiver photomixer due to the incoming THz wave, and this photocurrent then informs us of the amplitude of the electric field of the incoming wave. We can then also use a Fiber Stretcher attachment (which changes the optical path length difference between the transmitter and receiver) to get phase information in addition to the electric field amplitude.




Fig. 1: Transmission Setup: On the left we have the receiver photomixer and its collimating mirror, on the right we have the transmitter photomixer and its collimating mirror. Both photomixers are encased inside a silicon lenslet. The output from the transmitter is a beam that diverges by $12 - 15^\circ$. The collimating mirrors take the diverging beam and convert it into a plane wave (then take the plane wave and refocus the light into the receiver).

Laser Electronics

The figure below shows the laser electronics.

- HV amp: High voltage amplifier
- FS: Fiber stretcher
- DLC Smart: Digital Controller




Fig. 2: Laser electronics.

Laser Software

We control the laser through the TOPAS Software Graphical User Interface (GUI).




Fig. 3: Laser controller software.

Background: Reflection Measurement Setup

Figure 4 shows the reflection setup:




Fig. 4: Reflectometry setup.

See `Notebook2_Reflection.ipynb` for a walk through on carrying out reflection measurements.

Laser Spectroscopy: Measuring Materials' Optical Properties

Shreya Sutariya (shreyas@uchicago.edu), ERC 568 This is the hands-on component for the FLS lab.

Today's lab will center on reflection measurements of a common plastic called high-density polyethylene (HDPE). After measuring and plotting the measurement, we will plot its simulated performance and compare. Finally, we fit to get the exact index of refraction n of the material. (Time permitting we will also measure: LPE, MF-114, non-ARC material).

Today's Objectives:

- See `Notebook1_Laser_Spec_Background.ipynb` for material on understanding the measurement instrument and reflection setup.
- **Exercise 1:** Take reflection measurements of our sample.
- **Exercise 2:** Simulate what reflection of such a material should look like.
- **Exercise 3:** Fit our measurement to the model to get the index of refraction of our sample.
 - Bonus: We can understand a material more precisely by knowing its complex index of refraction, $\tilde{n} = n + ik$, where the imaginary component k determines the absorptivity of the material (and the real part n is what we normally refer to as the index of refraction). For code that can fit to both the real and imaginary parts, use the FLS Lab Computer to do additional fits using `Notebook3_Bonus_complex_index_fits.ipynb`.

Exercise 1: Reflection Measurements

As you take reflection measurements, think about the following:

- What can introduce noise in our measurement? What are ways we can reduce this?

Bonus:

- Think about the transmission setup from Figure 1 in `Notebook1`. What challenges and differences could come up in a transmission measurement?


Measurement Procedure:

In addition to measuring the sample's performance across a frequency range, we also have to take a calibration measurement to normalize the sample's measurement and get its reflectivity. We use a metal plate for this calibrating measurement, which serves as our 100% reflection.


For both the metal plate and the sample, we will get the electric field amplitude as it changes with respect to frequency. Then, to convert from amplitude to power (which scales as $|E|^2$), we square the sample and calibration plate amplitudes. To normalize our sample's measurement, we divide out the calibration plate: $R = (E_{\text{sample}}/E_{\text{plate}})^2$

Calibration measurement:


1. Set the sweep parameters in the TOPAS software:
 - Suggested values:
 - Start frequency: 60 GHz
 - End frequency: 180 GHz
 - Step size: 0.1 GHz



1. Set the current emitting frequency to 60 GHz.




1. Insert metal plate into sample holder; use an allen wrench and the screws to use washers to hold the plate in place.



2. Press 'Scan up' to start the sweep.



1. After the sweep is completed, choose the data columns to save and save your data in the `1ASchool_20220810` folder with the name `r_plate_1_60-180GHz.txt`.



1. Repeat these steps but now replacing the metal plate in step (3.) with the sample of choice (for the first measurement, use HDPE, for other measurements, ask me (Shreya) for more samples). Save the sample in the same folder as above with the sample's name, e.g. `r_hdpe_1_60-180GHz.txt`
2. To allow students to plot on their personal computers: Upload the folder and its files to the shared Google Drive folder: `cmb-s4_summer_school` using the link <https://drive.google.com/drive/folders/1dnvBpeqmDDmX2h6gwY8r5F7hnbjLUJ?usp=sharing>
 - In this case, download the files onto your personal computer and change the path later on to match this folder.

Exercise 2: Plotting our measurement

Here we plot the measurement from Exercise 1.

```
import numpy as np
import matplotlib.pyplot as plt
import scipy
import scipy.optimize as optimize
```

```
# group_folder_name = '1ASchool_20220810'
group_folder_name = '20220407'
name_sample = "alumina-1" # add a measurement number to keep track of things
```

```
# Next we have some useful functions/classes to smooth our data if needed, as well as
```

```
def smooth(x,window_len=11,window='hanning'):
    """smooth the data using a window with requested size.
    NOTE: length(output) != length(input), to correct this: return y[(window_len/2-1):-(window_len/2-1)]"""
    if x.ndim != 1:
        print('smooth only accepts 1 dimension arrays.')
    if x.size < window_len:
        print('Input vector needs to be bigger than window size.')
    if window_len<3:
        return x
    if not window in ['flat', 'hanning', 'hamming', 'bartlett', 'blackman']:
        print('Window is on of \'flat\', \'hanning\', \'hamming\', \'bartlett\', \'blackman\':')
    smp_r = [x[(window_len-1):0:-1],x[(2-(window_len-1)-1):]]
    if window == 'flat': #moving average
        wmp.ones(window_len,'d')
    else:
        w=eval('np.'+window+'(window_len)')
    ymp.convolve(w/w.sum(),s,mode='valid')
    return y

class fls_rt(object):
    def __init__(self, file_name):
        self.file_name = file_name
        return
    def load_data(self, path_ = "/Users/Shreya Sutariya/data/", col = 2): # col sees
        """
        Column indices correspond to: 0=freq set, 1=Thz photocurrent, 2=freq actual,
        """
        data = np.genfromtxt(path_+group_folder_name+'/'+self.file_name, skip_header=
        col_data = data[:,col]
        return col_data
    def get_rt(self, calibration_freq, calibration_ampitude, sample_freq, sample_amp):
        """check if the calibration arrays are the same lengths""" # check if the arrays
        len_calib = len(calibration_ampitude)
        len_sample = len(sample_ampitude)
        if len_calib != len_sample:
            # Interpolate for the calibration and sample arrays to match in length if
            f = scipy.interpolate.interpnd(calibration_freq, calibration_ampitude, b
            calibration_ampitude = f(sample_freq)
            rt = (sample_ampitude/calibration_ampitude)**2
            rt = (sample_ampitude/calibration_ampitude)**2
            return sample_freq, rt
    def get_rt_smooth(self, calibration_freq, calibration_ampitude, sample_freq, sam
        cal_f_smooth = smooth(calibration_freq, window_len = window_len_)
        cal_a_smooth = smooth(calibration_ampitude, window_len = window_len_)
        sam_f_smooth = smooth(sample_freq, window_len = window_len_)
        sam_a_smooth = smooth(sample_ampitude, window_len = window_len_)
        f, rt = self.get_rt(cal_f_smooth, cal_a_smooth, sam_f_smooth, sam_a_smooth)
        return f, rt
    def save_file(self, freq, refl, name_save, path_save = "/Users/Shreya Sutariya/dat
        header = "Frequency(GHz), Reflectivity"
        data = np.array([freq, refl])
        data = data.T
        outfile = open(path_save+name_save, 'w')
        outfile.write(header)
        np.savetxt(outfile, data)
        outfile.close()
        print("Your file is saved and stored in: ", path_save+group_folder_name+'/'
```

Fill in the file names.


```
### Insert file names here:
alumina1 = fls_rt("alumina2_60-120ghz.txt")
plate1 = fls_rt("plate2_60-120ghz.txt")
```

If doing this on your personal computer, use the last four lines with the `path_` argument and comment out the first four lines. Change the `path_` argument to match wherever you downloaded the files from Google Drive.

```
alumina_freq = alumina1.load_data(col=2) # frequency
alumina_ampitude = alumina1.load_data(col=3) # electric field amplitude
plate_freq = plate1.load_data(col=2)
plate_ampitude = plate1.load_data(col=3)
# alumina_freq = alumina1.load_data(path_ = "/Users/(Your Name)/Downloads/", col=2) #
# alumina_ampitude = alumina1.load_data(path_ = "/Users/(Your Name)/Downloads/", col=3)
# plate_freq = plate1.load_data(path_ = "/Users/(Your Name)/Downloads/", col=2) #
# plate_ampitude = plate1.load_data(path_ = "/Users/(Your Name)/Downloads/", col=3)
```

```
f, r = alumina1.get_rt(plate_freq, plate_ampitude, alumina_freq, alumina_ampitude)
f_smooth, r_smooth = alumina1.get_rt_smooth(plate_freq, plate_ampitude, alumina_freq,
alumina1.save_file(f_smooth, r_smooth, "example.txt", path_save = "/Users/Shreya Suta
Your file is saved and stored in: /Users/Shreya Sutariya/data/20220407/saved_files/ex
ample.txt
```

```
plt.figure()
plt.plot(f, r, alpha=0.5)
plt.plot(f_smooth, r_smooth)
plt.xlabel("Frequency [GHz]")
plt.ylabel("Reflectance")
plt.show()
```



Trouble Shooting: Standing Waves

You'll notice that there are quite a few wiggles at the peaks in the measurement.

- What do you think these could be due to? What we do know about standing waves?

A: Standing Waves:

- Where in the system could we imagine standing waves originating?
- What are methods to mitigate standing waves? (Hint: Look at materials on the optical bench)
- Why did we put Eccosorb (black absorbing foam) on the sample mount?

After the adjustment to the setup to try and mitigate standing waves, re-do the measurement steps from Exercise 1 (re-measure plate and sample, this time saving the files with the next measurement number, e.g. `plate2_60-180GHz.txt` and `hdpe2_60-180GHz.txt` if this is the second measurement of the day).

Then, replot the data using these two code cells:

```
# Insert file names here
hdpe2 = fls_rt("hdpe2_60-120ghz.txt") # change "sample_num" to sample name with the m
plate2 = fls_rt("plate2_60-120ghz.txt")
```

```
hdpe_freq2 = hdpe2.load_data(col=2) # frequency
hdpe_ampitude = hdpe2.load_data(col=3) # electric field amplitude
plate_freq2 = plate2.load_data(col=2)
plate_ampitude2 = plate2.load_data(col=3)
f2, r2 = hdpe2.get_rt(plate_freq2, plate_ampitude2, hdpe_freq2, hdpe_ampitude2) # u
f_smooth2, r_smooth2 = hdpe2.get_rt_smooth(plate_freq2, plate_ampitude2, hdpe_freq2,
plt.figure()
plt.plot(f2, r2, alpha=0.6, "C1")
plt.plot(f_smooth2, r_smooth2, "C1")
plt.plot(f, r, alpha=0.5, "C2")
plt.xlabel("Frequency [GHz]")
plt.ylabel("Reflectance")
plt.show()
```

Exercise 2: Simulating the Material Performance

We can model a slab of some material with parallel edges as a Fabry-Perot cavity. As waves travel in and reflect back from the second boundary, interference occurs with the incoming rays. Both rays have a difference in the optical path length traveled, which results in a phase difference, causing a periodic interference pattern.

We can determine what the exact spacing of this interference pattern should look like by knowing what the thickness d of the sample is, as well as the index of refraction of the material n , and the angle of incidence during our reflection measurement.

- Measure the material thickness in mm using the caliper provided.

- Determine the angle of incidence by (1) taking a bird's eye view photo of the setup and (2) using this site to measure the angle between the receiver and transmitter, going out to the sample. Divide this value by two.


Note: This code does not take into account any absorptive properties of the material!

```
def AiryR(params):
    """
    n0, n, freq, d, angleI, R = params
    """
    n0, n, freq, d, angleI, R = params
    R = ((n-1)/(n+1))**2
    c = 3e8
    lamb = c/freq
    omega = 2*np.pi*(c/lamb)
    k = n*omega/c
    k0 = n0*omega/c
    angleT = np.arcsin((n0/n)*np.sin(angleI))
    l = d*np.cos(angleT)
    s = 2*d*np.sin(angleT)/np.cos(angleT)
    a = s*np.sin(angleI)
    diff_phase = 2*k*l - k0*a
    T = (1 + 4*R/(1-R))**2*(np.sin(diff_phase/2))**2)**(-1)
    return 1-T
def functionAiryR(nl, constants, ff_):
    angleI_deg, thickness_mm = constants
    angleI_rad = np.radians(angleI_deg) # units: radians
    thickness_m = thickness_mm*1e-3
    # print("this is thickness:", thickness_m)
    ff_Hz = ff_*1e9 # units: Hz
    n0 = 1 # index of air: n0 | n1 | n0 where n1 is the index of the material
    params = [n0, n1, ff_Hz, thickness_m, angleI_rad]
    R = AiryR(params)
    return R
```

In the next cells, plot a simulated version of what interference pattern we expect to see:

```
# 1. Choose your frequency range:
freqs_sim = np.linspace(60, 120, 300) # start freq [GHz], end freq [GHz], number of p
```

```
# 2. Choose your parameters
# a. predicted material index n1
# b. angle of incidence, theta_i [deg]
# c. material thickness, d [mm]
plt.figure()
plt.plot(freqs_sim, functionAiryR(3.0, [13.5, 3.175], freqs_sim, '--', label='Sim 1')
# plt.plot(freqs_sim, functionAiryR(3.0, [13.5, 5], freqs_sim, '-', label='Sim 2') #
# plt.plot(freqs_sim, functionAiryR(3.0, [13.5, 10], freqs_sim, '-.', label='Sim 3') #
plt.plot(f_smooth, r_smooth, label='Measured')
plt.xlabel("Frequency [GHz]")
plt.ylabel("Reflectance")
plt.legend()
plt.show()
```



Uncomment out the extra simulation plotting lines and play around with different indices, thicknesses, etc.

- What do you think will happen to the interference pattern as the material thickness increases?
- How about as the index of refraction increases?

Plot our measurement data from Exercise 1 on top of the simulation and compare the two by eye. Play around with the simulation parameters to see how the two can match better.

Exercise 3: Fitting to the Index of Refraction

Using our measured data from Exercise 1 and the rough parameter estimates from Exercise 2, fit to the material's index of refraction.

```
def err(n1, constants, ff_, sampler):
    model = functionAiryR(n1, constants, ff_)
    return np.sum((samplerR - model)**2)
```

Enter the index guess below, as well as the angle of incidence and thickness.

Then, in the `X_fit = optimize.minimize(_` line, enter the frequency and reflection of the sample that you want to fit to.


```
p0 = 3. # Best n1 guess from Exercise 2
constants = [13.5, 0.125*25.4] # theta_i, incidence angle [deg], material thickness d [m]
# enter the measured frequencies and reflection data in the "arguments" (args):
X_fit = optimize.minimize(err, p0, args=(f_smooth, r_smooth), method='Nelder-M
```

```
n_fit = X_fit.x[0]
print("The result of fitting to the index of refraction is n1 = ", n_fit)
```

The result of fitting to the index of refraction is `n = 3.0319335937500007`

Now plot the simulation with this index of refraction and our measured data.

```
plt.figure()
plt.plot(f_smooth, functionAiryR(n_fit, [13.5, 3.175], f_smooth, '--', label='Fit')
plt.plot(f_smooth, r_smooth, label='Measured')
plt.xlabel("Frequency [GHz]")
plt.ylabel("Reflectance")
plt.legend()
plt.show()
```



Great! Now we have our fit.

- What could be done to improve the fit? (Remember, we haven't accounted for the imaginary component of the index of refraction here. How well do we trust our measurement of the material thickness? The angle of incidence?)


Success: You are done with `Notebook 2: Reflection`! You can now measure another material. If you finish early with the extra materials and have time, check out `Notebook 3: Bonus Complex Index fits`.

Additional materials:

- Low Pass Edge (LPE) filter from 180 to 250 GHz
- Alumina (aluminum oxide) disc from 60-120 GHz

Additional Material: Low Pass Edge (LPE) Filter

After repeating the measurement steps from Exercise 1 for the plate and LPE (take care with the LPE, it's more fragile), plot the results using the next two cells. The LPE is a metal-mesh filter, described in more detail here.



```
# Insert file names here
lpe = fls_rt("lpe3_180-250ghz.txt") # change "sample_num" to sample name with the mea
plate = fls_rt("plate3_180-250ghz.txt") # change plate number to match sample measure
```

```
lpe_freq = lpe.load_data(col=2) # frequency
lpe_ampitude = lpe.load_data(col=3) # electric field amplitude
plate_freq = plate.load_data(col=2)
plate_ampitude = plate.load_data(col=3)
f_lpe, r_lpe = lpe.get_rt(plate_freq, plate_ampitude, lpe_freq, lpe_ampitude) # un-
f_smooth_lpe, r_smooth_lpe = lpe.get_rt_smooth(plate_freq, plate_ampitude, lpe_freq,
plt.figure()
plt.title("LPE")
plt.plot(f_lpe, r_lpe, "C1", alpha=0.6, )
plt.plot(f_smooth_lpe, r_smooth_lpe, "C1")
plt.xlabel("Frequency [GHz]")
plt.ylabel("Reflectance")
plt.show()
```

```
OSError: Traceback (most recent call last)
<ipython-input-41-53015886f46b> in <module>
----> 1 lpe_freq = lpe.load_data(col=2) # frequency
      2 lpe_ampitude = lpe.load_data(col=3) # electric field amplitude
      3
      4 plate_freq = plate.load_data(col=2)
      5 plate_ampitude = plate.load_data(col=3)
<ipython-input-10-7e998e79653e> in load_data(self, path_, col)
      8 Column indices correspond to: 0=freq set, 1=Thz photocurrent, 2=freq a
      9 ctual, 3=amplitude, 4=phase
--> 10 data = np.genfromtxt(path_+group_folder_name+'/'+self.file_name, skip
header=1)
      11 col_data = data[:,col]
      12 return col_data
~\anaconda3\lib\site-packages\numpy\lib\newapi.py in genfromtxt(fname, dtype, comments,
delimiter, skip_header, skip_footer, converters, missing_values, filling_values, usec
ols, names, excludelist, deletechars, replace_space, autostrip, case_sensitive, default
fmt, unpack, usemask, loose, invalid_raise, max_rows, encoding)
   1747 fname = os.fspath(fname)
   1748 if isinstance(fname, str):
-> 1749 fid = np.lib._datasource.open(fname, 'rt', encoding=encoding)
   1750 else:
   1751     fid_ctx = contextlib.closing(fid)
~\anaconda3\lib\site-packages\numpy\lib\_datasource.py in open(path, mode, destpath, e
ncoding, newline)
   193
   194 ds = DataSource(destpath)
--> 195 return ds.open(path, mode, encoding=encoding, newline=newline)
   196
   197
~\anaconda3\lib\site-packages\numpy\lib\_datasource.py in open(self, path, mode, encod
ing, newline)
   533 encoding=encoding, newline=newline)
--> 535 else: raise IOError("%s not found." % path)
   536
   537
OSError: /Users/Shreya Sutariya/data/20220407/lpe3_180-250ghz.txt not found.
```

- Looking at this plot, high reflection means low transmission. Around which frequency does the LPE stop transmitting and instead reflect heavily?

References

On photomixers and THz lasers:

- [Info Sheet: Laser-based terahertz generation & applications](#), M. Lang, et al.
- [Review of resonating continuous wave terahertz systems and current application trends in terahertz domain](#), R. Safian, et al.

More on LPEs:

- [PowerPoint: Metal-mesh technology: a past and presetti view](#), L. Moncelis, et al.
- [A review of metal mesh filters](#), P. Ade, et al.